

Date of finalization:
July 23, 2024

@ Email:
info@emterm.com

Official website:
www.emterm.com

Corporate website:
www.watmedical.com

HeadTerm

Headache Migraine Relief Device

USER MANUAL

Headache Relief
Without Medication
Version NO.: US01

Return & Help Center
For any product failure not related to user error,
please contact us for a free replacement.

[Note]

- Relieve Migraine: 86% Effective in Reducing Migraine Pain Levels. The device is used to relieve migraine, not to eliminate them completely.
- Please clean the forehead and the gel side of the electrode before using the device.
- For the best experience, use the device while lying down and secure the electrode with the head strap if necessary.
- The device will automatically turn off when the battery is low. If it has not been used for an extended period, ensure it is fully charged before using it again.
- The device automatically stops treatment after 60 minutes of continuous use.
- Replace the electrodes as needed to provide optimal adhesive contact.

- For first-time use, it is recommended to start at a low gear level (1-2) to better adapt to the tingling sensation. Increase the gear as needed in subsequent sessions.
- For acute migraine treatment: Complete a 60-minute treatment session. If the migraine persists, start another 60-minute session.
- For migraine prevention: Accumulate a total of 60 minutes of use per day. Regular use over time yields better results.
- Ensure the device is turned off before removing it from the forehead.

Table of Contents
[QUICK START GUIDE] 01 [Risks and Benefits of Device Use] 21
[Name of Product] 02 [Warnings] 22
[Specification and Model] 02 [Precautions] 24
[Packaging Specification] 02 [Cleaning and Maintenance] 25
[Instructions] 03 [Disposal] 26
[Operation Guide] 11 [Operating Conditions] 26
[SPECIFICATIONS] 15 [Storage] 27
[Structure of Product] 16 [Guidance and Manufacturer's Declaration-Electrode] 18 [Electromagnetic Emission] 28
[Essential Performance] 18 [Labels and Marks] 29
[Indications for Use] 19 [EMC Information] 31
[Operation Principle] 19
[Date of Manufacturing and Shelf Life] 20
[Contraindications] 20

QUICK START GUIDE

HeadTerm

TENS device-HeadTerm 2
[Specification and Model]
VF-H12

[Packaging Specification]
The packaging specification for each model is one piece per case.
The product configuration list is shown in the following table:

Name of Accessory	Quantity
Main unit	1 set
Electrode	3 pads/per pack
User manual	1 set
USB-C cable	1 set
Head strap	1 set

Note: Electrode pack quantity varies by version.

Instructions:
1. Clean the forehead with a cleansing wipe, then dry it before using the device. This will help the electrode attach to the forehead securely.

2. Peel off the plastic protective layer from the gel side of the electrode.

3. With the gel side facing down, use a mirror for guidance and place the electrode one finger-width above the eyebrows.

4. The main unit magnetically attaches to the electrode.

5. Operating instructions
5.1. Turn on: Press and hold the control button for 2 seconds while the device is turned off.

5.2. Single click the control button for 1 gear increase.

5.3. Please adjust the gear to the level you feel most comfortable to achieve the best results.

5.4. Turn off: Press and hold the control button for 2 seconds while the device is turned on.

5.5. The head strap may be used to prevent the electrode from detaching.

6. Remove the device
6.1. Remove the main unit from electrode.
6.2. Peel the electrode sheet off the forehead and cover the gel with the plastic protection cover.

7. Clean the main unit with a cleansing wipe.

8. Store the main unit and electrodes in the storage box.

9. Please use CE or FCC certified adapter. It is recommended to use the computer USB port to charge the device.

Operation Guide:

Operation	Actions	Motor Feedback	Buzzer Feedback	Indicator Feedback
Turn on	Long press the control button for 2 seconds	Motor vibrates twice	/	Blinks periodically
Increase gear	Tap the control button for 1 gear intensity increase	/	1 x "beep" sound	Blinks periodically

Operation	Actions	Motor Feedback	Buzzer Feedback	Indicator Feedback
Decrease gear	Tap the control button twice	/	2 x "beep" sounds	Blinks periodically
Turn off	Long press the control button for 2 seconds	Motor vibrates twice	/	Light off
Battery level check (fully battery charged)	Click the control button in the turn off state	/	/	Remains on for 2 seconds
Battery level check (low battery level)	/	/	/	Blinks for 2 seconds

Operation Guide:

Status	Explanation	Motor Feedback	Buzzer Feedback	Indicator Feedback
Treatment completed	After 60 minutes of continuous use, the device automatically shuts down	Motor vibrates twice	/	Light off
In 9th gear	Maximum intensity, cannot increase any further	Motor vibrates twice	/	Blinks periodically
In 1st gear	Minimum intensity, cannot decrease any further	Motor vibrates twice	/	Blinks periodically

Status	Explanation	Motor Feedback	Buzzer Feedback	Indicator Feedback
Low battery warning	Activates every five minutes when under 20% battery	Motor vibrates three times	/	Blinks periodically
Charging	Charging the device with USB-C cable	/	/	Blinks periodically
Fully charged	/	/	/	Remains on
Unable to be turned on	Insufficient battery to be turned on	/	/	/

[Indications for Use]

The TENS device-HeadTerm 2 is indicated for the prophylactic treatment of episodic migraine in patients 18 years of age or older.

[Operation Principle]

By releasing the low frequency pulse with a particular frequency and reaching the advanced nerve center of cerebral cortex via the supraorbital nerve, device can adjust signals causing headache which come from the biological and physical channels to stop or postpone the transmission of headache signal to cerebral center. Meanwhile, the device can improve an individual's ability to cope with headache symptoms by activating the release of endorphins in the body.

[Structure of Product]
1. Structure of the product:
Main unit.

No.	Name	Function
1	Control button	Turns device on and off
2	Indicator	Decrease/increase intensity
3	Back (electrode)	Indicates operating status
3	Back (electrode)	Generates low frequency pulse

Figure 1. Device Illustration

[Risks and Benefits of Device Use]
Risks of device use:

- A stinging sensation may occur on the forehead.
- If you experience pain, try lowering the intensity by double-clicking the control button to decrease the level.

Benefits of device use:

- Prophylactic treatment of episodic migraine through non pharmacological means.
- No drug side effects.

Not adhering to the care regimen may result in ineffective treatment and potential injuries (see contraindications).

[Warnings]

- Do not use in the presence of electronic monitoring equipment (e.g., cardiac monitors, electrocardiogram (ECG) alarms).
- Do not use while bathing or showering.
- Do not use during any activity that can put you at risk of injury.
- EMC WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- EMC WARNING: Because of this device adjacent to or stacked with other devices should be avoided because it could result in improper operation. If such use is necessary, this device and the other devices should be observed to verify that they are operating normally.
- EMC WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device Connected.

[Cleaning and Maintenance]
HeadTerm 2 Main Device
Step 1: Turn off the device.
Step 2: Wipe the surface of the device with a cleaning wipe, especially the magnetic contact area connected to the electrode.
Step 3: Perform a visual inspection. If the surface is still contaminated, repeat step 2 until it is clean.
USB-C Charging Cable
Handle the cable carefully to avoid tangling. Keep the connectors at both ends of the cable clean. If they become dirty, clean them with a cleaning wipe.
HeadTerm 2 Electrodes
After use, reattach the protective film to the electrode pads to prevent dust and debris buildup. Store them in the original storage box to preserve adhesion and prolong their lifespan.

[Disposal]
Dispose device according to local state/county laws for electronic waste.
[Operating Conditions]
Normal operating conditions should comply with the following requirements:
Environment temperature: 5°C~+55°C (-4°F~+131°F)
Relative humidity: 10%~80%
Atmospheric pressure: 50~106 kPa
Indoor, dry and well-ventilated, free from corrosive substances
Do not put excessive pressure device or store under heavy items.

[Labels and Marks]

No.	Symbol	Meaning
1		Manufacturing date of product
2		Serial number of product
3		Note! Please check all documents attached
4		Information of manufacturer
5		BF type
6		Device contains a battery. Do not dispose in general waste stream/household trash.

[Guidance and Manufacturer's Declaration - electromagnetic immunity - for all EQUIPMENT and SYSTEMS]

The VF-H12 TENS device-HeadTerm 2 is intended for use in the electromagnetic environment specified below. The customer or the user of the VF-H12 TENS device-HeadTerm 2 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD)	IEC 61000-4-2	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the surface resistance should be at least 10 ⁹ Ω.
Electrostatic transient / burst	IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for differential mode	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage	IEC 61000-4-11	100% 0.5 Hz to 100 Hz 95% 100 Hz to 1 kHz 95% 1 kHz to 10 kHz 95% 10 kHz to 100 kHz 95% 100 kHz to 1 MHz 95% 1 MHz to 10 MHz 95% 10 MHz to 100 MHz 95% 100 MHz to 1 GHz	Main power quality should be that of a typical commercial or hospital environment. If the use of the VF-H12 TENS device-HeadTerm 2 requires constant operation during power mains interruption, it is recommended that the VF-H12 TENS device-HeadTerm 2 be powered from an uninterrupted power supply or a battery.
Power frequency magnetic field	IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: (1) Is the A.C. mains voltage prior to application of the test level.

[Date of Manufacturing and Shelf Life]
See the packaging.
Device shelf life: 36 months
[Contraindications]

- Presence of metallic or electronic implants in the brain, or a pacemaker.
- Acute inflammation, hemorrhagic tendency, arrhythmia, or epilepsy.
- Traumatic brain injury, head trauma, or maxillofacial surgery.
- Brain tumors, meningitis, or acute cerebrovascular stroke.
- You should check with your doctor if you do not understand or are unsure whether any of the above contraindications apply to you or your condition.

[Structure of Product]
1. Structure of the product:
Main unit.

No.	Name	Function
1	Control button	Turns device on and off
2	Indicator	Decrease/increase intensity
3	Back (electrode)	Indicates operating status
3	Back (electrode)	Generates low frequency pulse

Figure 1. Device Illustration

[Electrode]
The electrode acts as the interface between the main unit and the skin. The electrode has two contact zones that connect with the HeadTerm 2 through magnets.

Magnetic contact zone	Conductive gel

[Essential Performance]
Pulse repetition amplitudes: ±60 V(±3 V)
Pulse duration: 250 us (±20 us)
Pulse repetition frequencies: 60 Hz(±1 Hz)

[Appendix 1 - Guidance and manufacturer's declaration - electromagnetic immunity - for all EQUIPMENT and SYSTEMS]

No.	Guidance and manufacturer's declaration - electromagnetic immunity
1	The VF-H12 TENS device-HeadTerm 2 is intended for use in the electromagnetic environment specified below. The customer or the user of the VF-H12 TENS device-HeadTerm 2 should ensure that it is used in such an environment.
2	Immunity test
3	Compliance
4	Electromagnetic environment - guidance
5	RF emissions
6	RF emissions (CISPR11)
7	Class B
8	Harmonic emissions
9	IEC 61000-3-2
10	N/A
11	Voltage fluctuations / RF emissions
12	IEC 61000-3-3
13	N/A

[Appendix 2 - Guidance and manufacturer's declaration - electromagnetic immunity - for all EQUIPMENT and SYSTEMS]

No.	Guidance and manufacturer's declaration - electromagnetic immunity
1	The VF-H12 TENS device-HeadTerm 2 is intended for use in the electromagnetic environment specified below. The customer or the user of the VF-H12 TENS device-HeadTerm 2 should ensure that it is used in such an environment.
2	Immunity test
3	Compliance level
4	Electromagnetic environment - guidance
5	RF emissions
6	RF emissions (CISPR11)
7	Class B
8	Harmonic emissions
9	IEC 61000-3-2
10	N/A
11	Voltage fluctuations / RF emissions
12	IEC 61000-3-3
13	N/A

[Appendix 3 - Guidance and manufacturer's declaration - electromagnetic immunity - for EQUIPMENT and SYSTEM that are not LIFE SUPPORTING]

No.	Guidance and manufacturer's declaration - electromagnetic immunity
1	The VF-H12 TENS device-HeadTerm 2 is intended for use in the electromagnetic environment specified below. The customer or the user of the VF-H12 TENS device-HeadTerm 2 should ensure that it is used in such an environment.
2	Immunity test
3	Compliance level
4	Electromagnetic environment - guidance
5	RF emissions
6	RF emissions (CISPR11)
7	Class B
8	Harmonic emissions
9	IEC 61000-3-2
10	N/A
11	Voltage fluctuations / RF emissions
12	IEC 61000-3-3
13	N/A

[Appendix 4 - Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTEM that are not LIFE SUPPORTING]

Rated maximum output power	150 MHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.7 GHz
0.01	0.035	0.035	0.035
0.1	0.11	0.11	0.11
1	0.35	0.35	0.35
10	1.1	1.1	1.1
100	3.5	3.5	3.5

WAT Med
Designed in Canada.
WAT Medical Enterprise Ltd.
Add: Unit 170, 422 Richards St, Vancouver, BC, V6B 2Z4, Canada