


Wet-Stop³⁺

Select Vibrate (V)
Alarm (A), or Both (VA)



Insert batteries
correctly
#2

TO TURN OFF


1. Push the LED button.
2. Wipe sensor dry OR disconnect sensor cable.

TIPS

- Child should wear alarm every night until dry for at least 2-3 weeks.
- If alarm does not wake child, wake him/her.
- Keep a log of dry and wet nights (reward chart included).
- This is a team approach: stay positive!
- If alarm gets wet, open battery door and remove batteries. Let dry for several days. If this does not work then you will need a new alarm.
- Download and read The Complete Bedwetting Book at PottyMD.com or WetStop.com.

CAUTION

1. Not intended to be used as an infant training aid for potty training.
2. Not intended for children under four years of age, unless recommended by a physician.
3. Make sure that the cord does not become tangled around the user.
4. To prevent a choking hazard, do not place any components in the mouth at any time.
5. A soft rug or padding at the bedside will help prevent injury should your child roll out of bed upon awakening. Minimizing the elevation of the bed above the floor is advisable.
6. The Wet-Stop³⁺ alarm should only be used for one person.



CLEANING



Alarm is NOT WATERPROOF!

The alarm cannot be rinsed or submerged in water or urine.

Wipe alarm and sensor with dry cloth as needed.
Only sensor can be cleaned with soap and water.



One year warranty against manufacturer's defects in materials and workmanship and against breakage or failure under normal use. This warranty does not cover misuse or abuse.



1-865-584-6700



FREE DOWNLOAD

The Complete Bedwetting Book
wetstop.com



Instructions and rewards
wetstop.com



MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany

DOC. #WS103.2 English



CE Marking of Conformity



Type B equipment



See instructions for use



Battery