



Malem Enuresis Alarms Instruction Booklet

IMPORTANT
Please read instruction booklet
before use

Best Results are obtained when:

- 1. Child is motivated and wants to be dry.**
- 2. Child wakes up to the alarm when activated.**
- 3. Parents/Guardians help to wake up the child when the alarm sounds especially at the early stages of treatment.**
- 4. No lifting, no nappies and no withholding fluids.**

Alarms cure by speeding up the natural process

Contents

Setup for MO3, MO4/S, MO5 MO6 & MO17.....	3-5
Important/Hints	6
Batteries	7
Choice Selection	8
Sensors	8
Fault Finding	8
Guarantee	9
Repairs and Reconditioning	9
Cleaning	9
Electrical Safety and Disposal.....	9
Bed-Mat Enuresis Alarm - MO3, MO4/S, MO5, MO6 & MO17	10-12
Wander-Mat sensor (WM)	13
Dialysis-Mat sensor (DM)	13
All sensors.....	14
Model MO3 - Enuresis Alarm	15
Model MO4 - Ultimate 1 Enuresis Alarm.....	16
Model MO4S - Ultimate 1S Enuresis Alarm	17
Model MO5 - Ultimate 1 + Record Enuresis Alarm.....	18
Model MO6V - Vibrating Unit for Bedside Alarm.....	19
Model MO6 - Ultimate Bedside Alarm.....	20-21
Model MO17	22-23
Model MO8 - Sit Up Sensor	24
Standard, Easy-Clip™ and Bed-Mat sensors	25
Dialysis-Mat sensor (DM)	26
Wander-Mat sensor large and small (WM)	27
Wander-Mat inverter	28-29
Model PC9 - Queen's Square Bladder Stimulator.....	30
Model PC9-S - Queen's Square Bladder Stimulator (9 modes) ..	32


Malem Enuresis Alarms

Instructions for the MO3, MO4, MO5, MO6 & MO17 alarms. For Bed-Mat alarms see PAGES 10-12


1. Insert sensor plug into alarm socket making sure the **Easy-Clip™ Sensor** lever is lifted up (see diagram). Do not unplug unless for cleaning or replacement. (Remove battery isolation tag if present)
2. Carefully attach the alarm to pyjama jacket or night dress close to the shoulder and thread cord inside top as shown. For small children, loop and tie cord to shorten length. Attach to material before wearing garment if preferable.
3. For Vibrating alarms, attach on the inside of close fitting top or T-shirt close to the collar bone.
4. For VCA (Light Emitting), attach in a position on client's or child's clothing where it can readily be seen by carer or nursing staff. Alternatively, attach alarm to the end of the bed where it can be seen by nursing staff and secure sensor/mat under bed sheet.
5. For **Easy-Clip™ Sensor**: Lift sensor clip lever, attach to the outside of close fitting underwear, at the location most likely to get wet first, by inserting material inside sensor jaws at the front and closing lever. (Dispose of insulating tag if present)



Lowering the lever without any dry material in the jaws will trigger the alarm and prevent the alarm from resetting. Lift the lever and press the reset button.




6. For **Standard Sensor**: Cut off one end of the absorbent pad (not supplied) and insert the sensor in between the layers. Anchor with the crocodile clip to prevent sensor from being pulled out. Uncover the backing strip and stick the pad inside a pair of close fitting pants, insuring the sensor is located where it will come in contact with the flow of urine.



**PLEASE NOTE THAT
SENSOR & CLIP ARE
GOLD PLATED.**

Boys who prefer not to use pads can sandwich the sensor between **TWO PAIRS** of cotton briefs (not boxer shorts) or place within fold of Y-Front, and anchor by locating the crocodile clip on the outside below the waistband of the pants. Adhesive tape may also be used to fasten the sensor plate correctly in the area of urine flow. The sensor plate and crocodile clip **must not** be in direct contact with the skin. Avoid nylon/polyester underwear.

7. For **Bed-Mat Sensor**: See instruction on page 10.



Please make sure the sensor plate and the crocodile clip are attached between the two pairs of pants and do not come in contact with the skin.

8. Upon waking up to the alarm, try to stop any more urine from being released from the bladder. Remove the Standard sensor and wipe dry with a cloth. Remove the Easy-Clip™ sensor and keep the lever open. Do **NOT** unplug the cable from the body of the alarm. Stop the sound or vibration by pressing the small reset button, use the lavatory to completely empty the bladder, replace the sensor into a fresh pad or clip to a dry pair of pants as before. **Always test alarm by wetting sensor or closing the Easy-Clip™ lever before use.**
Parental help in waking up children when the alarm sounds is reassuring and helpful especially during the first few nights.

Vibrating ‘PCT’ - Personal Continence Trainer

- a. The vibration signals the start of urination - retrain the bladder by tensioning the pelvic floor muscles. With practice, it is possible to stop urination until the toilet is reached.
- b. Note down the times of urination as they occur. Gradually, over a period of time, a pattern will be apparent and the bladder can be emptied before wetting.

Light Emitting ‘VCA’ - Visual Continence Alarm

- a. Can be used to give visual indication of urination so that the client is changed to prevent prolonged exposure of skin to urine with associated risk and complications.
- b. Can be used to develop a toileting regime by tabulating the client’s urination times. A pattern can be established and the need for a bedpan or use of toilet anticipated.
This can result in a significant reduction in management cost as well as an improvement in client’s moral.

MO4 & MO4S Ultimate 1 Alarm

The choice of Sound/s Only, Sound/s and Vibration or Vibration Only is selectable by sliding the lever inside the battery compartment to setting 1, 2, or 3 accordingly. See diagram on page 16 and 17.

MO5 Ultimate 1 + Record Alarm

For choice selection of sound, vibration or recording see diagram on page 18. For MO17 see page 22 and 23.

Important

- **Rapid wakening on hearing or feeling the alarm is vital for success. Most children, including non-bedwetters, tend to sleep soundly and heavily at home. Better results are obtained if Parents/Guardians help to wake the child up on hearing the alarm. The child must wake up sufficiently to be able to recall the event in the morning.**
- Use the progress chart. Great encouragement is derived by seeing steady improvement.
- Always praise the user for a dry night; never scold or chastise for a wet night. Encouragement is the key to success. You can not help a child who does not want to help him or her self. Try and find out why they do not want to be helped, it could be some other problem.
- To avoid relapse use the alarm for a minimum of two weeks after dryness has been achieved. At the same time, increase fluid intake before sleep to achieve over learning and avoid possibility of relapsing at a later stage. Avoid fizzy drinks, tea, coffee & milk.

Hints

- **The Reset Button will stop the alarm from sounding only if the Standard Sensor plate is dry or the Easy-Clip™ lever is LIFTED OPEN and REMOVED from the underwear.**
- **Do not unplug the sensor from the Alarm.**
- There is no on/off switch that could allow the alarm to be switched off and ignored. Current drain on the batteries is minute except when the alarm is sounding.
- The sensor is positively locked to the alarm to prevent accidental disconnection. However, the sensor can be unlocked and detached for cleaning or replacement. **Do not unplug the sensor from the alarm during normal use.**
- It is intentional that one of the sounds of the Eight Tone Alarm (8T) does not stop immediately upon pressing the reset button.
- **REGULARLY CLEAN the Easy-Clip™ sensor in soapy water, rinse and shake dry to prevent “fluff” and urine accumulation inside the jaws.**

Batteries (remove isolating tag if present)

Batteries are provided free of charge by ourselves during manufacture and units are tested prior to packaging. Batteries are safely secured within the alarm. **They should not be accessed by children.**

Batteries should be replaced only by adults.

For model MO3 button size **Alkaline batteries** AG13 or equivalent. L1154, SP675H, V357, LR44, H1154, MR44, RW52, PX675, 675RP, HG13, SR44 or SG13. See diagram on page 15.

PLEASE DO NOT USE ZINC AIR BATTERIES (also known as hearing aid batteries) AS THEY CONTINUOUSLY DISCHARGE REGARDLESS OF BUZZING AND ARE UNABLE TO SUSTAIN POWER TO THE ALARM. (They have a sticker when new and tiny holes on the +ve side).

Do not allow ribbon to slip between batteries.

For **Ultimate 1 Alarm, Model MO4 & MO4S**, use 2 x AAA Alkaline batteries (see diagram on page 16 & 17). For **Ultimate 1 + Record Alarm, Model MO5 and MO17** use 3 x AAA Alkaline batteries (see diagram on page 18 & 23). For Model **MO6**, use 4 x AA Alkaline batteries. To change the batteries, remove the battery cover. Take special care to observe correct polarity when replacing all the batteries (see diagram on page 20). Make sure batteries are not loose within the compartment. Do NOT mix old and new batteries.

Batteries do tend to deteriorate during use and storage and it is advisable to test the unit before use; this can be simply done by moistening the finger and touching the sensor or closing the Easy-Clip™ lever. The alarm should always sound, vibrate or light up. If the unit does not respond, replace the batteries. In use, the batteries should not be expected to last more than 1-2 months, so it is advisable to test the unit regularly. Maximum continuous operating time is about 60 minutes per set of button batteries and up to 10 hours for the AAA Alkaline batteries.

Important - please remove batteries and clean sensor if the alarm is not going to be used for some time.

Choice Selection

To change setting combination for MO4 & MO5 remove battery cover and slide lever inside battery compartment (see diagram on page 16 and 18). Position 1 = Sound/s Only. Position 2 = Sound/s + Vibration. Position 3 = Vibration Only. For MO4S see page 17.

To select message or sound for MO5 & MO6 (see diagram on pages 18, 20 & 21). To record a new message make sure the selection switch is in message position, then press record button and speak into the microphone at the back of the alarm. Maximum recording time is 10 seconds. **Recording a message from Mum or telephone tones is very effective. For MO17 see page 22.**

Sensors

Sensors are an important component and must always be clean and dry. The crocodile clip is gold plated and should be used to anchor the sensor in the correct position. **Wash clean with soapy water and rinse. Remove any material from jaws of Easy-Clip™ sensor. Dry before use. See page 14.**

Sensor Wire

The wire that connects the sensor to the alarm unit is a vitally important link. The wire is sufficiently flexible to be comfortable in use but sufficiently robust to ensure long life. It will not however tolerate abuse, being bitten or chewed. **Do not use the wire to pull the sensor out of the pad or underwear - always undo the clip.** Do not clip the crocodile clip to the sensor wire. Do not unplug the sensor from the body worn alarms.

Fault Finding

- a) **Unit fails to sound:-** replace batteries, ensure correct polarity. Always use “fresh batteries”. Do not use Zinc/Air batteries (the ones with the small holes in the positive side and sticker when new). Make sure ejector ribbon has not slipped between batteries. Make sure battery contacts are clean and free from deposits. Make sure batteries are not loose and are making good contact.
- b) **Unit will not turn off:-** sensor plate is wet (Dry sensor). Sensor is faulty or has a damaged wire replace sensor. Easy-Clip™ sensor lever is closed or material is wet, lift lever, remove material & press the reset button. Easy-Clip™ is dirty and requires cleaning see page 3. Batteries are nearly exhausted please replace batteries. Alarm is or has been flooded.
- c) **Intermittent buzzing:-** faulty sensor. Loose or flat batteries in the battery compartment. Damaged wire.

Guarantee

Under normal usage the Malem Alarm Unit is guaranteed for a period of **three months** from purchase. **This excludes batteries, sensors, Wander-Mats, Dialysis-Mats, Bet-Mats, flooding, misuse or accidental damage.** If the alarm should malfunction during this period please contact us **before sending**. Please remove the batteries and decontaminate the unit/sensor before returning. (see below). Do not return flooded or contaminated alarms.

Repairs and Reconditioning

If your alarm requires repair, we will exchange it for a new one (repairs are not undertaken for hygiene reasons). Please contact us by telephone, e-mail or in writing, for the cost of this service. Please decontaminate the unit before return.



Cleaning

Alarm Unit: Wipe clean alarm and safety pin using mild detergent or skin disinfectant. Use a damp cloth and dry thoroughly. **NEVER IMMERSER IN WATER.**


Standard sensor and all Mats: Rinse with detergent and dry thoroughly.

Easy-Clip™ sensor: **Lift lever and regularly clean in soapy water, rinse and shake dry. If the alarm is to be reused by another client/patient/child then use a new sensor or Bed-Mat. The National Guidelines states one child, one sensor.** To replace the sensor, unlock the sensor by lifting the small tag and gently pulling the plug out of the alarm.

Electrical Safety

Designed to meet all Department of Health and British Standard BS5724 Part 1, RE1004/03 and IEC performance, construction and safety requirements. Maximum current across the sensor is less than 10µA. BSI certificate No.: 115627. BSI EMC test certificate No.: EZ000505. Electrical Safety Symbol  Type BF.  Marked.

Disposal

At the end of the alarm working life, please dispose of the product at a special local Authority Differentiated Waste Collection Centre. 

Malem Bed-Mat Enuresis Alarm

Instructions for MO3, MO4, MO5, MO6 and MO17


- Under no circumstances should any of the alarms be attached to the users' clothing when using a Bed-Mat.
- Place a clean and dry Bed-Mat correct side up over the mattress in the area which is normally wetted and gently tie down.
- Cover the Bed-Mat with a Clean Dry **Cotton sheet and/or Bath Towel** and make the bed up as usual. Avoid nylon/polyester sheets as they can cause 'false' alarms due to increased perspiration. Adjust cover thickness over the mat to suit.
- Do not place plastic sheets directly under the Bed-Mat as condensation will cause the alarm to trigger. (See diagram page 11)
- Attach the Bed-Mat wire to the chosen Malem Alarm and place the alarm beside the bed apart from the Vibrating alarm which can be placed under the pillow.
- Please tuck in the Bed-Mat wire so that it is not in the way for getting in and out of bed.
- Upon waking to the alarm, try to stop any more urine from being released and use the lavatory to completely empty the bladder.
- The alarm sound can be stopped **ONLY** when the Bed-Mat wire is unlocked from the alarm and the reset button is pressed.
- The Bed-Mat must be wiped clean with detergent and dried then replaced in a freshly made bed, as before, for reuse.
- Adjust the loudness of the sound with the volume control. (MO6 only)
- To test the alarm, place a damp cloth across the foil tracks on the Bed-Mat.
- Bed-Mat sensitivity can be reduced by sliding switch on side of MO6 Alarm. (See diagram on page 20)

Under NO circumstances should any of the alarms be attached to the users' clothing when using a Bed-Mat.




Holes for tying Bed-Mat to the mattress.

If using a plastic protective mattress sheet place adhesive tape over the Bed-Mat lower contacts.



Use a BATH TOWEL to cover the Bed-Mat and cotton sheet. Adjust towel thickness to suit.



Pillow

Bedside alarm or other small Malem Alarm.

Place cotton sheet over Bed-Mat.

Tie Bed-Mat around mattress.

Tuck wire from Bed-Mat to alarm away neatly to prevent tripping or entanglement.

Bedside table

Important

- Rapid waking on hearing the alarm is vital for success but parents believe that bed-wetting children are very heavy sleepers so better results are obtained if parents help to wake the child when the alarm sounds especially during the first few nights.
- To avoid relapse use the alarm for at least 2 weeks after dryness has been achieved. At the same time, increase fluid intake before sleep to achieve over learning and avoid possibility of relapsing at a later stage. Avoid fizzy drinks, tea, coffee and milk.
- Use the progress chart. Great encouragement is derived by seeing steady improvement.
- Always praise the user for a dry night; never scold or chastise for a wet night. Encouragement is the key to success.

Hints

- The reset button will stop the alarm from sounding only if the Bed-Mat is dry or disconnected.
- The Bed-Mat positively locks to the alarm to prevent accidental disconnection. However, it must be unlocked before attempting to reset the alarm.
- Please use a **thicker cotton sheet or place a bath towel** above the Bed-Mat if sweating is triggering the alarm prematurely. Avoid plastic and man-made fibres as they may cause static discharge and false alarm. Reduce alarm sensitivity. (See diagram on page 20)
- Treat the Bed-Mat with care. Avoid folding it or jumping over it.
- Please replace the Bed-Mat with a new one if it has deteriorated and use a new one with each child. (National Guidelines)
- Please note that sleeping directly on an uncovered Bed-Mat or touching it by hand may trigger the alarm.
- Do not place plastic sheets directly under the Bed-Mat as condensation may trigger the alarm.

If used with Wander-Mat Sensor (WM):

- a) **Activation by applied pressure:** place the Wander-Mat (WM) in the appropriate position depending on the way you wish to use it. The Wander-Mat can be placed beside the bed or near the door so that treading on the Mat will activate the alarm. Please place a door mat over the Wander-Mat so as not to trip over the Wander-Mat. Plug into the alarm and place away safely. See page 27.
- b) **Activation by removal of pressure:** place the Wander-Mat (WM) over the mattress but under the sheet, in a position so that the user is laying over it or place on a chair so that the user is sitting on it. Plug the Wander-Mat wire into the 'inverter' then plug the inverter wire into the alarm and place away safely. If the user gets out of bed or out of the chair, the alarm will be activated. See page 27.
- c) Once the alarm is triggered, the LED will be flashing. The buzzer sound can be silenced by pressing the reset button on the top of the alarm. If the sensor remains wet or active then the alarm will sound again within 5 seconds.

If used with Dialysis-Mat Sensor (DM):

- a) Select the correct side of the Mat with shiny tracks up.
- b) Place the Mat over the chair arm.
- c) Plug the Mat wire into the alarm, tuck the wire and place the alarm out of the way in a safe place. See page 26.
- d) Test the Mat by placing a slightly wet tissue over the tracks. The alarm should sound.
- e) If sweating is likely to activate the alarm and sound the alarm cover the Mat surface with tissue paper. Adjust thickness to suit.
- f) Rest patients arm over the Mat.
- g) Wipe Mat clean and dry before use again.

Sensors

Easy-Clip™: Length 79 cm. Keep lever up when not in use. Clean sensor regularly. See page 3 and 25.

Standard: Length 79 cm. Keep clean and wipe dry. The sensor plate and crocodile clip MUST NOT be in direct contact with the skin. See page 4 and 25.

Bed-Mat: 54 x 42 cm. Wire length 220 cm. Keep clean and wipe dry. Always cover with cotton sheet as well as a bath towel. Under no circumstances should the alarm be attached to user's clothing when used with a Bed-Mat. See page 25.

Wander-Mat (Large): 75 x 40 cm. Wire 150 cm. Use pressure on or pressure off (for pressure off use with 'inverter'). See page 27, 28 and 29.

Wander-Mat (Small): 40 x 37 cm. Wire 150 cm. Use pressure on or pressure off (for pressure off use with 'inverter'). See page 27, 28 and 29.

Dialysis-Mat: 22 x 17 cm. Wire 100 cm. Keep clean and wipe dry. See page 26.


Malem Enuresis Alarm

MODEL MO3

Lift up lever to attach alarm to garment


Slide folded garment material between the body of the alarm and lever grip

Lower lever down to lock garment to the alarm



Carefully attach to garment near shoulder

Reset button



Twist slotted screwdriver to slide battery cover

Make sure battery ejection ribbon does NOT slide between batteries

Please observe correct polarity when replacing batteries

Batteries are AG13 or equivalent. L1154, SP675H, V357, LR44, H1154, MR44, RW52, PX675, 675RP, HG13, SR44 or SG13.

Use Alkaline batteries ONLY
No Zinc/Air please

Malem Ultimate 1 Enuresis Alarm

MODEL MO4


Lift up lever to attach alarm to garment

Slide folded garment material between the body of the alarm and lever grip

Lower lever down to lock garment to the alarm



Reset button



Carefully attach to garment near shoulder

Twist slotted screwdriver to slide battery cover

Use 2 x AAA Alkaline batteries

Please observe correct polarity when replacing batteries

Sensor locking socket

Combination Selection Switch

Slide to:

- 1 For Sound/s only
- 2 For Sound/s and vibration
- 3 For Vibration only

Use Alkaline batteries ONLY

Malem Ultimate 1S Enuresis Alarm

MODEL MO4S

Selection Switch

Slide to:

1 = Single Sound*

8 = 8 Sounds**

*See sound selection switch for choice.

**Sound selection switch must be on position 1.




Lower lever down to lock garment to the alarm

Lift up lever to attach alarm to garment



Slide folded garment material between the body of the alarm and lever grip

Reset button




Combination Selection Switch

Slide to:

S For Sound/s only

V&S For Sound/s and vibration

V For Vibration only



Carefully attach to garment near shoulder

Sensor locking socket

For 8 sound choice

Only switch position 1 must be on (rest in off position)

Single sound selection switch use a very small screwdriver or paper clip to activate switch - (be careful)

Twist slotted screwdriver to slide battery cover


Use 2 x AAA Alkaline batteries
Please observe correct polarity when replacing batteries



Use Alkaline batteries ONLY

Malem Ultimate 1 + Record Enuresis Alarm


MODEL MO5



Use Alkaline batteries ONLY

Vibrating Unit for Malem Ultimate Bed-side Alarm


MODEL MO6V



Use Alkaline batteries ONLY

Malem Ultimate Bed-side Alarm

MODEL MO6




Use Alkaline batteries ONLY

★ Single sound selection switch must be on sound 1 (Top switch) for 8 sounds selection to be active


Combination Selection Switch

Slide to: 1 = 8 Sounds* 1 = Message
 2 = Single Sound 2 = Sounds



Single sound selection switch
 On ←→ Off
 Select only one

- ★ Single sound
- 1
 - 2
 - 3
 - 4
 - 5
 - 6
 - 7
 - 8




Test Button

This LED is on when recording a new message as long as the record button is pressed

Combination Selection Switch

Slide to:
 1 = LED + Sound
 2 = LED only
 3 = No LED or Sound

Press to record (Max 10 secs recording time).
 To record a message the message/sound Combination Selection Switch must be in position 1 (message)




Malem™ Ultimate Multi-Choice Alarm

MODEL MO17

1. **SELECTION** - Press repeatedly to select one of 3 messages, one of 20 melodies or one of 8 sounds (See Selection Switches)
2. **STOP** - Alarm will reactivate after 5 sec if sensor remains active
3. **RESET** - Stop sound or LED (As long as sensor is deactivated)

Vibration Selection Switch



SELECTION SWITCHES


Recording

Make sure selection switches are on 'message' and 'sequence'. Press and hold the record button found under the battery cover for 3 seconds. Wait for red LED to light up, then speak clearly into the microphone at the back of the alarm. Maximum recording time is 10 seconds.

Message Melody Sound LED


Random Lock Sequence

The lock position will fix the selected choice



In-line recording socket

LED Blinks when alarm is triggered and when recording a message




Speaker

Dimensions 73mm x 69mm x 26mm Weight 95g including battery

Use Alkaline batteries ONLY

Malem™ Ultimate Multi-Choice Alarm


MODEL MO17



Use Alkaline batteries ONLY

Malem Sit Up Sensor

MODEL MO8



The Alarm will activate if the users movement make the sensor stand upright.

For best results:

1. Attach alarm inside a Pyjama Top pocket or
2. Attach on the inside of a T-Shirt Collar or
3. Clip alarm to an elasticated bandage around the thigh.

To Stop Alarm:


1. Place Sit Up Sensor in horizontal position and press the reset button on the side.
2. Remove Sit Up Sensor from alarm and press the Reset button.



Slide folded garment material between the body of the alarm and lever grip

Lift up lever to attach alarm to garment

Lower lever down to lock garment to the alarm




Carefully attach to garment near shoulder

Push Sensor up into locking socket to secure


Sit Up Sensor can be used with all Malem Alarms

Activation of the alarm using different sensors


1. Activation by Easy-Clip™ Sensor or Standard Sensor



2. Activation using a Bed-Mat:




3. Activation using a Dialysis-Mat




Wander-Mat with inverter (Large and Small)

4. Activation by applying pressure:




5. Activation by removing pressure:

(user **must be** lying in bed or sitting over WM in a chair)



Please note that pressure must be applied to the Mat BEFORE attempting to set up the connection.

Malem Wander-Mat Inverter front:




Please plug in the inverter wire between the Wander-Mat and the transmitter.
Battery test: press the button, LED will glow. If LED does not glow then change the battery. (Remove battery isolator if present)

The inverter should be used **ONLY** if you wish to use the Wander-Mat to detect **removal** of pressure.

Please do not immerse inverter or Mat in water.

Use Alkaline batteries ONLY

Malem Wander-Mat Invertor back:



Please observe correct polarity. Contacts clean and battery not loose inside compartment. **Pull metal contact to make good connection.** Use only 23A 12V Alkaline battery.

Use Alkaline batteries ONLY

Malem Queen's Square Bladder Stimulator

MODEL PC9

To activate
press pad against skin or
use switch

Push button down for ON
Push again for OFF

Pull Down

Use 2 x AA
Alkaline batteries
Please observe
correct polarity
when replacing
batteries

How to use

1. Place the
Stimulator on
the lower
abdomen
(above the
pubic bone)
before commencing
passing urine and
apply gentle pressure
to activate the vibration.
(If you are unable to
exert sufficient pressure
to activate, then use the on/off
button on the back)

2. Keep the stimulator vibrating while urine is being passed
and for about one to two minutes after you feel you have finished.

Model PC9



Adjust rope
length by tying
a knot

Use Alkaline batteries ONLY

Malem Queen's Square Bladder Stimulator

MODEL PC9-S (9 modes)

To activate
press pad against skin or
use switch

Pull Down

Use 3 x AA
Alkaline batteries
Please observe
correct polarity
when replacing
batteries

How to use

1. Place the Stimulator on the lower abdomen (above the pubic bone) before commencing passing urine and apply gentle pressure to activate the vibration. (If you are unable to exert sufficient pressure to activate, then use the on/off button on the back)

2. Keep the stimulator vibrating while urine is being passed and for one to two minutes after you feel you have finished. Please select different mode for most effective result. Press & hold mode button for two seconds (while stimulator is on) to activate. Button remains active for sixty seconds.

Push button down for ON
Push again for OFF

Mode
Button

MODE
Press & hold
for 2 seconds

OFF
ON

Model PC9-S
CE

Adjust rope
length by tying
a knot

Use Alkaline batteries ONLY

Disclaimer

Malem™ Alarms are designed as a helpful and versatile aid. It is not intended to replace human care or supervision and the manufacturer, Malem Medical, is not responsible for any consequences resulting from the use, misuse or failure of this device.

The manufacturer reserves the right to change specifications without prior notice.



Alarms cure by speeding up the natural process

Malem Medical, 10 Willow Holt, Lowdham
Nottingham NG14 7EJ, England, UK

Tel: + 44 (0)115 966 4440 Fax: + 44 (0)115 966 4672

E-mail: malem@malem.co.uk

Web: www.malem.co.uk Web: www.malemmedical.co.uk



© Malem Medical 2012