

"I can't be awake
around the clock ..."



**VALUE
PACK***

Babysense 2



"... so I use
my Babysense."

Sleep with greater peace of mind
thanks to Oricom's Babysense2.

**Redesigned compact styling with the
same proven technology and operation
as the original Babysense2.**

Babysense2 monitors the breathing movement rate of your sleeping baby. Should these breathing movements fall below a certain rate or not be detected, Babysense2 will alert you audibly.

Oricom's Babysense2 – the next level of care.

- Included as a medical device on the Australian Register of Therapeutic Goods (ARTG No. 97479).
- Babysense2 continually monitors the breathing movement rate of your baby.
- The control unit connects to sensor pads placed under your baby's mattress.
- An audible alarm is activated if your baby's breathing movement slows to a rate of less than 10 breaths per minute, or cannot be detected for more than 20 seconds.

 **oricom®**

Connecting you now.

www.oricom.com.au

*BS2SC710 Includes Babysense2 Infant Breathing Movement Monitor and Secure710 Video Monitor.

Babysense2 is intended for use as an Infant Apnoea Alarm.
ALWAYS READ THE LABEL AND USER GUIDE, USE ONLY AS DIRECTED.